

GECKO

Version 1.4.0 - Detailed Release Notes

Main Menu Updates

- **UI Consistency Improvements:** All scrollable panels now include scrollbars. Buttons are grouped and clearly labeled to improve discoverability and overall visual appeal.
- **Channel Search Bar:** A search field allows users to quickly filter channels when inspecting or selecting export variables. Due to space limitations, it's currently only a single search bar that will simultaneously filter the channels in the inspect channel area, and also the export channels area.
- **Session State (session_state.json):** The application now saves the user's working state inside the .kinarm export subfolder. When reopening the program, users can resume where they left off, including selected trial, export channels, and selected event markers. Note that it will overwrite the currently selected channels, so if you have Gaze_X and Gaze_Y selected, and you load a new .kinarm file you were working on previously, the Gaze_X and Gaze_Y might be removed if you didn't have them selected in the previous session.
- **Clear Cache Button:** Previously labeled 'Clear Interpolation Cache'. Users can now choose between clearing interpolation data, labeling order preferences, or the saved session state. Clearing any will re-prompt the window, when necessary, like it was their first time using the program.
- **Help Button:** Added under Utilities. Contains usage guidance, feature explanations, disclaimers, and file path information for exported CSV/JSON files.

Labeler GUI Updates

- **Plot Stability Fixes:** Resolved visual glitch where a smaller duplicate plot would flash or remain on screen. Fixed the legend positioning to eliminate excessive flickering during labeling.
- **Maximized Windows:** Labeler and interpolation windows now open maximized for better visibility and precision during manual annotation.
- **Centered Controls:** Buttons should be centered at the bottom of the interface for cleaner layout. This may or may not change depending on screen resolution.

- **Trial Progress Indicator:** The top banner now displays progress in the format 'current_trial / total_trials' to provide orientation during labelling.
- **Custom Labeling Order:** During the first time labeling, users can define the order in which fixation, pursuit, and saccade labeling occurs. This preference is persistent until cleared.
- **Edit Other Button:** On the summary screen--allows marking general or undefined gaze events. These export as 0 in the gaze_event column and appear as light orange in the plot.
- **Y-Axis Label Clarification:** Updated label to 'Distance (cm)' and necessary units to cm.
- **Erase Mode:** Toggleable erase functionality allows users to remove labeled segments by selecting a start and end frame. Provides controlled correction without restarting the trial or relying on the undo button.
- **±20 Frame Buffer:** Users may now select slightly outside trial bounds for precision. Selections are automatically clipped to valid frame ranges (e.g., -10 to 50 becomes 0 to 50).